



REF SI-11142

9 Fr. distal lumen | 2 Lumen | 10 cm device length | .035 inch dia. spring-wire guide

# MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 8 Fr. Catheters

**Contents:**

- 1: Two-Lumen Central Venous Access Device<sup>1</sup>: 10 cm Radiopaque Polyurethane with 9 Fr. Distal Lumen, Integral Hemostasis Valve, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 20.7 cm
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape
- 1: Dust Cap: Vented
- 3: Dust Cap: Non-Vented
- 1: Stopcock: 4-way High-Flow
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

<sup>1</sup>Licensed under US Patent No. 5,250,038.

- fr** Jeu d'accès veineux central à deux lumières MAC avec valve hémostatique intégrée, à utiliser avec les cathéters de 7 - 8 Fr.
- de** Einführbesteck für zweilumigen, MAC-zentralen Venenkatheter mit Hämostaseventil mit Kathetern der Größen 7 - 8 Fr.
- it** Set di accesso venoso centrale a due lumi MAC con valvola emostatica incorporata per l'uso con cateteri da 7 - 8 Fr.
- pl** Dwukanałowy zestaw MAC do uzyskiwania dostępu do żył centralnych, z integralnym zaworem hemostatycznym, do użycia z cewnikami o rozmiarze 7 - 8 F
- pt** Conjunto para Acesso Venoso Central MAC com Lúmen Duplo, Válvula de Hemostase Integral, para utilizar com Cateters de 7 - 8 Fr.
- ru** Комплект для катетеризации главных вен MAC с двухпросветным катетером и несъемным гемостазным клапаном, для использования с катетерами 7 - 8 Fr.
- sl** Komplet za centralni venski dostop MAC z dvema svetlinama z vgrajenim hemostatskim ventilom za uporabo z katetri 7 – 8 Fr.
- es** Conjunto para Acceso Venoso Central MAC con Lúmen Duplo, Válvula de Hemostase Integral, para utilizar con Cateters de 7 - 8 Fr.
- sv** MAC tvåkanaligt set för centralt ventillträde med integrerad hemostasventil för användning med 7 - 8 Fr. katetrar
- tr** MAC İki Lümenli Santral Venöz Erişim Seti, Entegre Hemostaz Valfi, 7 - 8 Fr. Kateterler ile kullanılmak üzere


Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

⚠ California Prop. 65  
Warning: Cancer and Reproductive Harm  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.

Lumen	Priming Volume*		Flow Rate† (mL/hr)
	(mL)		
Distal (9 Fr.) no catheter	1.65		33000
Distal (9 Fr.) with 8 Fr. catheter	0.90		10500
Proximal (12 Ga.) no catheter	0.60		13000

\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

EC REP Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
2797 Co. Westmeath, Ireland

**Teleflex®**  
www.teleflex.com  
II-11142-131C (2021-04)



REF SI-11142

9 Fr. distal lumen | 2 Lumen | 10 cm device length | .035 inch dia. spring-wire guide

MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 8 Fr. Catheters

