

REF AJ-01671

LATEX-FREE

PICC Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

ARROW®

Special Order Product

Contents:

- One: Indwelling Catheter: 16 Ga. x 27-1/2" (70 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Injection Site Cap, Integral Suture Wing, Contamination Guard
- One: Spring-Wire Guide: .035" (.89 mm) dia. x 39-1/2" (100 cm) PTFE-Coated (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™¹
- One: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- One: Fastener: Catheter Clamp
- One: Pressure Transduction Probe
- One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow Raulerson Spring-Wire Introduction Syringe²
- One: Clamp: Catheter
- One: Tissue Dilator
- One: Scalpel: #11
- One: Stopcock: 3-way
- One: UserGard® Needle-Free Injection Hub
- One: Injection Site Cap

¹ US Patent No. 5,484,419.

² Licensed under US Patent Nos. 4,813,938 and 5,045,065.

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

This product contains **NO** natural rubber latex.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

1101671-101E (11/16)

REF
AJ-01671

16^{Ga.} 70 cm
catheter
length

.035 inch dia.
spring-wire
guide

Peripherally Inserted Central Catheter Set with Blue FlexTip®

ARROW
INTERNATIONAL
Special Order Product

REF
AJ-01671

16^{Ga.} 70 cm
catheter
length

.035 inch dia.
spring-wire
guide

Peripherally Inserted Central Catheter Set with Blue FlexTip®

Peel to
Open
LATEX-FREE