



PICC Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

REF **AJ-01671**

16^{Ga.} | **70** cm catheter length | **.035** inch dia. spring-wire guide

Contents:

- | | |
|---|---|
| 1: Indwelling Catheter: 16 Ga. x 70 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Integral Suture Wing, Contamination Guard | 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm |
| 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 39-1/2" (100 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 2: Dust Cap: Non-Vented |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: SecondSite™ Adjustable Hub: Fastener |
| 1: Pressure Transduction Probe | 1: SecondSite™ Adjustable Hub: Catheter Clamp |
| 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Scalpel: #11 |
| | 1: Stopcock: 3-way |



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.70	1300

* Priming volumes are approximate and are done without accessories.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Rx only. **Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Manufacturer:
Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

teleflex®
www.teleflex.com

IJ-01671-101G (2021-04)



**PICC Peripherally Inserted Central Catheter Set
with Blue FlexTip® Catheter**

REF **AJ-01671**

16^{Ga.} | **70** cm catheter length | **.035** inch dia. spring-wire guide