



Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

REF **CK-03000-B**

7 Fr. | **2** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 3 mL Luer-Lock Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle


All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.

- fr** Kit de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Kit per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- pt** Kit de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- es** Juego de cateterización venosa central de dos luces con catéter Blue FlexTip
- sv** Uppsättning för tvåkanalig central venkatetrering med Blue FlexTip-kateter



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.40	3800
Proximal (16 Ga.)	0.38	3600

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA



IK-03000-135A (01/19)



Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

REF **CK-03000-B**

7 Fr. | **2** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide