



Central Venous Catheterization Kit

REF AK-04225

16^{Ga.} | 20^{cm} catheter length | .025^{inch dia.} spring-wire guide

Contents:

- 1: Indwelling Catheter: 16 Ga. x 20 cm Radiopaque Polyurethane with Removable Inner Plastic Cannula
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Local Anesthetic Component
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep^{®1} One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: Adapter: Arrow-Johans™ RAECG
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway[®] Disposal Cup
- 1: Drape: 24" x 40" (61 cm x 102 cm) with 2-5/8" (6.7 cm) fenestration, with adhesive
- 1: Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.



Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.


California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.26	4200

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



DuPont™ Tyvek[®] packaging

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex[®]
www.teleflex.com

IK-04225-100L (12/17)



Central Venous Catheterization Kit

REF AK-04225

16^{Ga.} | 20^{cm} catheter length | .025^{inch dia.} spring-wire guide