

Product No. AK-04400

Latex-Free

Central Venous Catheterization Kit

ARROW®

Catheter Diameter: 16 Ga., 1.7 mm

Priming Volume: 0.32 cc

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 cc.

Flow Rate: 3200 cc/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Contents:

- One: Indwelling Catheter: 16 Ga. x 12" (30 cm) Radiopaque Polyurethane with Side Holes, Integral Suture Wing
- One: Spring-Wire Guide: .032" (.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End — "J" Tip on Other)
- One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- One: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 cc Luer-Slip Syringe
- One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 cc Luer-Slip Syringe
- One: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 cc Luer-Lock Syringe
- One: Tissue Dilator
- One: Adapter: Arrow-Johans™ RAECG
- One: Clamp: Catheter
- One: Fastener: Catheter Clamp
- One: CSR Wrap
- One: Drape: 24" x 36" with 4" fenestration
- One: Scalpel: #11
- Three: Gauze Pads: 2" x 2"
- Four: Gauze Pads: 4" x 4"
- Two: Prep Sponge Swabs
- One: Suture: 3-0 Silk with Straight Needle
- One: Tray: Prep


Rx only.

Fluid path components are non-pyrogenic. Contains no medication.

This product contains **NO** natural rubber latex.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA



IK-04400-1001 (3/07)



Product No.
AK-04400

16 Ga.

30 cm catheter length

.032 inch dia. spring-wire guide

Central Venous Catheterization Kit

Peel to Open
LATEX-FREE