

REF AK-04660-S

# Central Venous Catheterization Kit

# ARROW®

Catheter Diameter: 24 Ga., 0.8 mm

Priming Volume: 0.10 mL

Priming volumes are approximate and are done without the injection cap.

Flow Rate: 690 mL/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Contents:

One: . . . . . Indwelling Catheter: 24 Ga. x 5" (12 cm) Radiopaque Polyurethane  
One: . . . . . Spring-Wire Guide: .018" (.46 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End – "J" Tip on Other)  
One: . . . . . Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW  
Two: . . . . . Syringes: 3 mL Luer-Lock  
One: . . . . . Clamp: Catheter  
One: . . . . . CSR Wrap  
One: . . . . . Drape: 18" x 26" with 3" fenestration  
One: . . . . . Scalpel: #11  
Three: . . . . . Gauze Pads: 2" x 2"  
One: . . . . . Prep Sponge Swab

Rx only.

**Warning:** Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

This product contains **NO** natural rubber latex.

Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 US

TECFLEX Teleflex Medical  
IDA Business and Technology Park  
Athlone, Ireland

www.arrowintl.com



- CZ** Centrální žilní katetrizační kit
- DK** Kateteriseringskit til centrale vener
- NL** Centraal-veneuze katheterisatiekit
- F** Kit de cathétérisme veineux central
- D** Besteck für zentralen Venenkatheter
- GR** Kit κεντρικού φλεβικού καθετηριασμού
- H** Centrális vénás katéterező készlet
- I** Kit per cateterismo della vena centrale
- PL** Pełny zestaw do cewnikowania żył centralnych
- P** Kit de Cateterização Venosa Central
- RUS** Комплект для катетеризации главных вен
- SK** Súprava pre centrálnu venóznú katetrizáciu
- E** Juego de cateterización venosa central
- S** Uppsättning för central venkatetrisering

IK04660131A (10/10)

**ARROW**  
INTERNATIONAL

REF  
AK-04660-S

24<sup>Ga.</sup> 12<sup>cm</sup>  
catheter  
length

.018<sup>inch dia.</sup>  
spring-wire  
guide

Central Venous Catheterization Kit

Peel to  
Open