

Product No. AK-04700

LATEX-FREE

Central Venous Catheterization Kit

ARROW®

Catheter Diameter: 14 Ga., 2.2 mm

Priming Volume: 0.35 cc

Priming volumes are approximate and are done without the injection cap. Injection priming volume is 0.17 cc.

Flow Rate: 6900 cc/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Contents:

One: Indwelling Catheter: 14 Ga. x 8" (20 cm)
Radiopaque Polyurethane with
Side Holes, Integral Suture Wing

One: Spring-Wire Guide: .032" (.81 mm) dia.
x 17-13/16" (45 cm) (Straight Soft Tip
on One End — "J" Tip on Other)

One: Catheter: 18 Ga. x 2-1/2" (6.35 cm)
Radiopaque over 20 Ga. RW Introducer Needle

One: Fastener: Catheter Clamp

One: Introducer Needle: 18 Ga. x 2-1/2"
(6.35 cm) XTW

One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)

One: Injection Needle: 25 Ga. x 1" (2.54 cm)

Two: Syringe: 5 cc Luer-Slip

One: Clamp: Catheter

One: Vessel Dilator

One: Drape: 18" x 26" with 3" fenestration

One: CSR Wrap

Four: Gauze Pads: 4" x 4"

Three: Gauze Pads: 2" x 2"

One: Scalpel: #11

Two: Prep Sponge Swabs

One: Suture: 3-0 Silk with Straight Needle

One: Tray: Prep

One: SharpsAway® Disposal Cup


Caution: USA Federal law restricts this device to sale by or on the order of a physician.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Contains no medication. Single use — do not resterilize. Sterilized by ethylene oxide.

This product contains **NO** natural rubber latex.

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA



IK-04700-100H, Rev. 1 (6/11)

ARROW®
INTERNATIONAL

Product No.
AK-04700

14^{Ga.}

20^{cm}
catheter
length

.032^{inch dia.}
spring-wire
guide

Central Venous Catheterization Kit

Peel to
Open
LATEX-FREE