

REF AK-09801

ARROW®



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheter

8.5 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other)
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Syringe: 3 cc Luer-Lock
- 2: Syringes: 5 cc Luer-Slip
- 1: Suture: 3-0 Braided Silk with Straight Needle
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Drape: 24" x 36" with 4" fenestration
- 1: CSR Wrap
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 1: Obturator: 8 Fr.
- 1: Scalpel: #11
- 2: Prep Sponge Swabs
- 1: Tubing: 20" Extension
- 1: Tray: Prep
- 1: SharpsAway® Disposal Cup

- fr** Kit d'introduction de gaines par voie percutanée avec valve hémostatique/port latéral intégrés pour emploi avec des cathéters 7 - 7,5 Fr.
- de** Perkutanes Schleuseneinführbesteck mit eingebauter Hämostaseventil/ Seitenanschluß-Baugruppe zum Gebrauch mit Kathetern der Größen 7 - 7,5 Fr.
- it** Kit per introduzione di guaina percutanea con valvola emostatica/bocchetta laterale incorporate per l'uso con cateteri da 7 - 7,5 Fr.
- pt** Kit de Introdução de Bainha Percutânea com Válvula de Hemostase/Porta Lateral, para utilizar com Cateteres 7 - 7,5 Fr.
- es** Juego intoductor de vaina percutánea con válvula hemostática y orificio lateral integral para uso con catéters de 7 a 7,5 Fr.
- sv** Upsättning för perkutan skyddsintroducerare med integral hemostasventil/sidoöppning för bruk med katetrar mellan 7 och 7,5 Fr.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.

Caution: This Product Contains Dry Natural Rubber.



STERILE EO



CE REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
2797

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

IK-09801-130B (2021-04)

REF AK-09801

ARROW®



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheter

8.5 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide