

REF REF AK-09801-A

ARROW



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheter

8.5 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other)
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Syringe: 3 cc Luer-Lock
- 2: Syringes: 5 cc Luer-Slip
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Drape: 24" x 36" with 4" fenestration
- 1: CSR Wrap
- 7: Gauze Pads: 4" x 4"
- 1: Obturator
- 1: Scalpel: #11
- 2: Prep Sponge Swabs
- 1: Suture: 3-0 Silk with Straight Needle
- 1: Tray: Prep
- 1: SharpsAway® Disposal Cup

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

This product contains dry natural rubber.

- fr** Kit d'introduction de gaines par voie percutanée avec valve hémostatique/port latéral intégrés pour emploi avec des cathéters 7 - 7,5 Fr.
- de** Perkutanes Schleuseneinführbesteck mit eingebauter Hämostaseventil/ Seitenanschluß-Baugruppe zum Gebrauch mit Kathetern der Größen 7 - 7,5 Fr.
- it** Kit per introduzione di guaina percutanea con valvola emostatica/bocchetta laterale incorporate per l'uso con cateteri da 7 - 7,5 Fr.
- pt** Kit de Introdução de Bainha Percutânea com Válvula de Hemostase/Porta Lateral, para utilizar com Cateteres 7 - 7,5 Fr.
- es** Juego intoductor de vaina percutánea con válvula hemostática y orificio lateral integral para uso con catéteres de 7 a 7,5 Fr.
- sv** Upsättning för perkutan skyddsintroducerare med integral hemostasventil/sidoöppning för bruk med katetrar mellan 7 och 7,5 Fr.



STERILE EO



CE 2797
 EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
 www.teleflex.com

IK-09801-132B (2021-04)

REF REF AK-09801-A

ARROW



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheter

8.5 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide