



Pediatric Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

REF **AK-14412**

4 Fr. | **2** Lumen | **13** cm catheter length | **.018** inch dia. spring-wire guide

Contents:

- | | |
|---|--|
| <ul style="list-style-type: none"> 1: Two-Lumen Indwelling Catheter: 4 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm | <ul style="list-style-type: none"> 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm 2: Dust Cap: Non-Vented 1: SecondSite™ Adjustable Hub: Catheter Clamp 1: SecondSite™ Adjustable Hub: Fastener 1: SharpsAway® Disposal Cup 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration 1: Scalpel: #11 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) 1: Suture: 3-0 Braided Silk with Straight Needle |
|---|--|

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Sterile. Fluid path components are non-pyrogenic.



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight




Keep dry



Do not use if package is damaged



Not made with natural rubber latex



Catheter Cross Section

Lumen	Priming Volume *(mL)	Flow Rate† (mL/hr)
Distal (22 Ga.)	0.25	550
Proximal (22 Ga.)	0.31	425

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

DuPont™ Tyvek®
sterile

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IK-14412-130F (5/17)



Pediatric Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

REF **AK-14412**

4 Fr. | **2** Lumen | **13** cm catheter length | **.018** inch dia. spring-wire guide