



Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

REF AK-15402

4 Fr. | **2** Lumen | **8** cm catheter length | **.018** inch dia. spring-wire guide

Contents:

- | | |
|--|---|
| <ul style="list-style-type: none"> 1: Two-Lumen Indwelling Catheter: 4 Fr. x 8 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW 1: Injection Needle: 25 Ga. x 1" (2.54 cm) 1: 3 mL Luer-Lock Syringe 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe 1: Tissue Dilator 1: Tissue Dilator 1: 5 mL 1% Lidocaine HCl Solution 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution with Hi-Lite Orange™ Tint | <ul style="list-style-type: none"> 2: Dust Cap: Non-Vented 1: SecondSite™ Adjustable Hub: Catheter Clamp 1: SecondSite™ Adjustable Hub: Fastener 1: SharpsAway® Disposal Cup 1: Drape: 18" x 26" with 3" fenestration 1: Local Anesthetic Components 1: Scalpel: #11 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) 1: Suture: 3-0 Braided Silk with Straight Needle <p>¹A registered trademark of CareFusion or one of its subsidiaries.</p> |
|--|---|


DuPont™ Tyvek® HDPE packaging



Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

All contents sterile unless otherwise indicated. Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr.)
Distal (22 Ga.)	0.25	540
Proximal (22 Ga.)	0.30	950

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®

www.teleflex.com

IK-15402-130C (1/17)



Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

REF AK-15402

4 Fr. | **2** Lumen | **8** cm catheter length | **.018** inch dia. spring-wire guide