

**REF AK-25703-A**

**ARROW**



**Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Arrowg+ard Blue® Catheter**

**7** Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

**Contents:**

- 1: Arrowg+ard Blue® Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.

**Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**



**Lumen**

Distal (16 Ga.)  
Medial (18 Ga.)  
Proximal (18 Ga.)

**Priming Volume\* (mL)**

0.40  
0.40  
0.40

**Gravity Flow Rate† (mL/hr)**

2720  
1260  
1390

\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



STERILE EO



Not made with natural rubber latex



Store between 20-25°C (68-77°F).

**REF AK-25703-A**

**ARROW**



**Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Arrowg+ard Blue® Catheter**

**7** Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Manufacturer:  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)

IK-25703-115A (2020/04)