



HANDS-OFF® Two-Lumen Central Venous Catheter with Blue FlexTip® for use only with Arrow One-Piece Percutaneous Sheath Introducer System

REF **HO-14702**

7 Fr. | **2** Lumen | **30** cm catheter length | **23** cm insertable catheter length with Cath-Gard®

Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 30 cm, Cath-Gard® Catheter Contamination Shield
- 2: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.


- fr** Cathéter veineux central HANDS-OFF à deux lumières avec Blue FlexTip uniquement pour emploi avec le système d'introduction de gaine par voie percutanée monobloc Arrow
- de** Zweilumiger zentraler HANDS-OFF-Venenkatheter mit Blue FlexTip ausschließlich zur Verwendung mit dem einteiligen Einführsystem für perkutane Schleusen von Arrow
- it** Catetere per vena centrale HANDS-OFF a due lumi con Blue FlexTip per l'uso esclusivo con il sistema di introduzione di guaina percutanea monopezzo Arrow
- pl** Dwukanałowy cewnik HANDS-OFF do cewnikowania żył centralnych z końcówką Blue FlexTip, do użycia tylko z jednoczęściowym systemem do przezskórne wprowadzania koszulki Arrow
- pt** Cateter Venoso Central HANDS-OFF de Lúmen Duplo com Blue FlexTip, para utilizar apenas com o Sistema de Introdução de Bainhas Percutâneas Arrow de Peça Única
- ru** Двухпросветный катетер для главных вен HANDS-OFF с Blue FlexTip, для применения исключительно с системой однокомпонентной чрескожной капсулы-интубатора Arrow
- sl** Centralni venski kateter HANDS-OFF z dvema svetlinama z Blue FlexTip je samo za uporabo z uvajalnim sistemom z enodelnim perkutanim tulcem Arrow
- es** Catéter para cateterismo venoso central de dos luces HANDS-OFF con punta Blue FlexTip sólo para uso con el sistema introductor de vaina percutánea Arrow de una pieza
- sv** HANDS-OFF tvåkanalig central venkateter med Blue FlexTip, endast för användning med Arrows odelade perkutana hylsinförarsystem
- tr** HANDS-OFF İki Lümenli Santral Venöz Kateter, Blue FlexTip ile, sadece Arrow Tek Parçalı Perkütan Kılıf İntroduser Sistemi ile kullanılmak içindir



Rx only.

Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.



Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA



Teleflex®
 www.teleflex.com

10-14702-135A (2020-12)



HANDS-OFF® Two-Lumen Central Venous Catheter with Blue FlexTip® for use only with Arrow One-Piece Percutaneous Sheath Introducer System

REF **HO-14702**

7 Fr. | **2** Lumen | **30** cm catheter length | **23** cm insertable catheter length with Cath-Gard®