



PICC Two-Lumen Peripherally Inserted Central Catheter with Blue FlexTip® Catheter Seldinger Access Set

REF **PR-05042**

4 Fr. | **2** Lumen | **50** cm catheter length | **.018** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 4 Fr. x 50 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: GlideThru™ Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure
- 1: Dressing: STATLOCK® Catheter Stabilization Device

†A registered trademark of C. R. Bard, Inc.

All components are CE 0086 unless otherwise noted.




- fr** Jeu d'accès Seldinger pour cathéter central inséré périphériquement CCIP à deux lumières avec cathéter Blue FlexTip
- de** PEZK-Besteck für zweilumigen, peripher eingeführten, zentralen Katheter für Seldinger-Zugang mit Blue FlexTip
- it** Set di accesso Seldinger al catetere centrale a due lumi a inserimento periferico con catetere Blue FlexTip
- pt** Conjunto de Cateter Central de Introdução Periférica PICC por Técnica de Seldinger, com Lúmen Duplo com cateter Blue FlexTip
- es** Conjunto de acceso Seldinger para catéter de vena central de inserción periférica (PICC) de dos luces con catéter Blue FlexTip
- sv** Set för PICC tvåkanalig perifert införd central kateter med Seldinger-tillträde med Blue FlexTip-kateter

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (22 Ga.)	0.30	130	740
Proximal (22 Ga.)	0.40	140	1010

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.
 ‡ Pump flow rates determined at maximum pump pressure of 10 psig.



Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA



IR-05042-131A (04/19)



PICC Two-Lumen Peripherally Inserted Central Catheter with Blue FlexTip® Catheter Seldinger Access Set

REF **PR-05042**

4 Fr. | **2** Lumen | **50** cm catheter length | **.018** inch dia. spring-wire guide