



## PICC Two-Lumen Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter Seldinger Access Set

REF **PR-05052**

**5** Fr. | **2** Lumen | **50** cm catheter length | **.018** inch dia. spring-wire guide

**Contents:**

- 1: Two-Lumen Indwelling Catheter: 5 Fr. x 50 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: GlideThru™ Peel-Away Sheath: 5 Fr. x 4" (10 cm) Radiopaque over 5 Fr. Dilator
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 1: Skin Protectant Prep Pad (non-sterile solution)
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure
- 1: Dressing: STATLOCK® Catheter Stabilization Device

- fr** Jeu d'accès Seldinger pour cathéter central inséré périphériquement CCIP à deux lumières
- de** PEZK-Besteck für zweilumigen, peripher eingeführten, zentralen Katheter für Seldinger-Zugang
- it** Set di accesso Seldinger al catetere centrale a due lumi a inserimento periferico
- pt** Cateter Central de Introdução Periférica PICC por Técnica de Seldinger, com Lúmen Duplo
- es** Conjunto de acceso Seldinger para catéter de vena central de inserción periférica (PICC) de dos luces
- sv** Set för PICC tvåkanalig perifert införd central kateter med Seldinger-tillträde

<sup>1</sup>A registered trademark of C. R. Bard, Inc.


All components are CE 0086 unless otherwise noted.



**Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.35	575	2450
Proximal (20 Ga.)	0.32	160	1025

\* Priming volumes are approximate and are done without accessories.  
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.  
 ‡ Pump flow rates determined at maximum pump pressure of 10 psig.

**Manufacturer:**  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA



**Teleflex®**  
 www.teleflex.com

IR-05052-140A (01/19)



## PICC Two-Lumen Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter Seldinger Access Set

REF **PR-05052**

**5** Fr. | **2** Lumen | **50** cm catheter length | **.018** inch dia. spring-wire guide