

Contents:

- One: Sheath: 9 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
 - One: Spring-Wire Guide: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other)
 - One: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
 - One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 cc Luer-Slip Syringe
 - One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 cc Luer-Slip Syringe
 - One: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 cc Luer-Lock Syringe
 - One: Suture: 3-0 Braided Silk with Straight Needle
 - One: 5 ml HCl, 1% Lidocaine Solution
 - One: 1 Applicator Pouch, 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution
 - One: Dust Cap: Non-Vented
 - One: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
 - One: Drape: 24" x 36" with 4" fenestration
 - One: Local Anesthetic Components
 - Two: Gauze Pads: 2" x 2"
 - Five: Gauze Pads: 4" x 4"
 - One: Obturator: 8 Fr.
 - One: Scalpel: #11
 - One: Stopcock: 4-way High-Flow
 - One: Tubing: 20" Extension*
 - One: SharpsAway® Disposal Cup
- ¹ A registered trademark of Medi-Flex Hospital Products Corporation.

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

Not made with natural rubber latex.

±Precaution: Contains Phthalate: DEHP.

* Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

INFECTION PROTECTION FEATURE:

- **Chloraprep®*** One-Step skin antiseptic, 2% chlorhexidine gluconate in 70% isopropyl alcohol



Product No. NR-09903-S

Percutaneous Sheath Introducer Kit

with Integral Hemostasis Valve/Side Port for use with 7.5 - 8 Fr. Catheters



Product No. **NR-09903-S** **9**^{Fr.} **10** cm sheath length **.035** inch dia. spring-wire guide
Percutaneous Sheath Introducer Kit

Product No. **NR-09903-S** **9**^{Fr.} **10** cm sheath length **.035** inch dia. spring-wire guide
Percutaneous Sheath Introducer Kit

NR-09903-100K (3/16)

Peel to Open