



# ARROW® PICC Set

## REF PR-34041-HPHNM

### Contents:

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>1: Indwelling TaperFree™ Catheter: 4 Fr. x 40 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamp, and Contamination Guard</li> <li>1: GlideThru™ Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator</li> <li>1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)</li> <li>1: Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil (CE 0120)</li> <li>1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW</li> <li>1: Syringe: 10 mL Luer-Lock</li> <li>1: Skin Protectant Prep Pad</li> <li>1: Dust Cap: Non-Vented</li> </ul> | <ul style="list-style-type: none"> <li>1: SecondSite™ Adjustable Hub: Catheter Clamp</li> <li>1: SecondSite™ Adjustable Hub: Fastener</li> <li>1: Catheter Trimmer</li> <li>1: Safety Scalpel: #11</li> <li>1: Flow Rate Information Card with Injection Log</li> <li>1: Patient ID Card</li> <li>1: Chart Sticker</li> <li>1: Patient Information Booklet</li> <li>1: Checklist/CLIP Sheet</li> <li>1: Paper Tape Measure</li> <li>1: Dressing: STATLOCK®<sup>1</sup> Catheter Stabilization Device</li> </ul> |
|--|---|

<sup>1</sup>A registered trademark of C. R. Bard, Inc.  
All components are CE 0086 unless otherwise noted.



Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland  
**CE**  
0086

**Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com