



ARROW® PICC Set

REF PR-35041-HPHNM

Contents:

- | | |
|--|---|
| 1: Indwelling TaperFree™ Catheter: 4 Fr. x 50 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamp, and Contamination Guard | 1: SecondSite™ Adjustable Hub: Catheter Clamp |
| 1: GlideThru™ Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator | 1: SecondSite™ Adjustable Hub: Fastener |
| 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends) | 1: Catheter Trimmer |
| 1: Guidewire, Marked: .018" (0.46 mm) dia. x 80 cm Hydrophilic-Coated Nitinol with Soft Tip Tungsten Coil (CE 0120) | 1: Safety Scalpel: #11 |
| 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW | 1: Flow Rate Information Card with Injection Log |
| 1: Syringe: 10 mL Luer-Lock | 1: Patient ID Card |
| 1: Skin Protectant Prep Pad | 1: Chart Sticker |
| 1: Dust Cap: Non-Vented | 1: Patient Information Booklet |
| | 1: Checklist/CLIP Sheet |
| | 1: Paper Tape Measure |
| | 1: Dressing: STATLOCK® ¹ Catheter Stabilization Device |

¹A registered trademark of C. R. Bard, Inc.
All components are CE 0086 unless otherwise noted.



Caution



Consult
instructions for use



Do not reuse



Do not resterilize



Sterilized
by ethylene oxide



Keep away
from sunlight



Keep dry



Do not use
if package is damaged



Not made with
natural rubber latex

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
CE
0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com