



Arterial Catheterization Kit

REF **SAC-00520-CA**

20^{Ga.} | **5** cm catheter length | **.021** inch dia. spring-wire guide

Contents:

- 1: Indwelling Catheter: 20 Ga. x 5 cm Radiopaque PEBA
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
- 1: Syringe: 3 mL Luer-Lock
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Needle Holder

- 1: Safety Scalpel: #11
- 1: Checklist Sheet
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®² 10 cm x 12 cm
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of 3M Company.



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA
Teleflex®
www.teleflex.com

IS-00520-103A (09/18)



Arterial Catheterization Kit

REF **SAC-00520-CA**

20^{Ga.} | **5** cm catheter length | **.021** inch dia. spring-wire guide