

REF PS-01651

ARROW



PICC Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

Contents:

- 1: Catheter: 16 Ga. x 50 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Peel-Away Catheter: 14 Ga. x 2-1/2" (6.35 cm) Radiopaque over 15 Ga. TW Introducer Needle and 3 mL Luer-Slip Syringe
- 1: Skin Protectant Prep Pad (non-sterile solution)
- 1: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Dressing: STATLOCK® Catheter Stabilization Device

*A registered trademark of C. R. Bard, Inc.

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

cs	Set s periférně zaváděným centrálním katetrem (PICC) s katetrem Blue FlexTip
da	PICC perifert indført centralt katetersæt med Blue FlexTip kateter
nl	PICC Perifeer in te brengen centrale katheterset met Blue FlexTip katheter
fr	Jeu de cathéter central inséré périphériquement CCIP avec cathéter Blue FlexTip
de	PEZK-Besteck für peripher eingeführten, zentralen Katheter mit Blue-FlexTip-Katheter
el	Σετ περιφερικά εισαγόμενου κεντρικού καθετήρα PICC με καθετήρα Blue FlexTip
hu	PICC periferiálisán bevezethető, centrális katéterező felszerelés Blue FlexTip katéterrel
it	Set per catetere centrale a inserimento periferico con catetere Blue FlexTip
pl	Wprowadzany obwodowo zestaw typu PICC do cewnikowania naczyń centralnych z cewnikiem z końcówką Blue FlexTip
pt	Conjunto de Cateter Central de Introdução Periférica PICC com Cateter Blue FlexTip
ru	PICC набор центрального катетера для периферийного введения с катетером Blue FlexTip
sk	Sada periférne zavádzaného centrálného katétra (PICC) s katétrom Blue FlexTip
es	Conjunto de catéter de vena central de inserción periférica (PICC) con catéter Blue FlexTip
sv	Set för PICC perifert införd central kateter med Blue FlexTip-kateter

Diameter: 16 Ga., 1.7 mm	Priming Volume* (mL)	Flow Rate† (mL/hr)
Lumen		
Distal (16 Ga.)	0.64	1850

* Priming volumes are approximate and are done without accessories.

† Flow rates are with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IS-01651-130D (09/19)

REF PS-01651

ARROW



PICC Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

Contents:

- 1: Catheter: 16 Ga. x 50 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Peel-Away Catheter: 14 Ga. x 2-1/2" (6.35 cm) Radiopaque over 15 Ga. TW Introducer Needle and 3 mL Luer-Slip Syringe
- 1: Skin Protectant Prep Pad (non-sterile solution)
- 1: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Dressing: STATLOCK® Catheter Stabilization Device

*A registered trademark of C. R. Bard, Inc.

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

cs	Set s periférně zaváděným centrálním katetrem (PICC) s katetrem Blue FlexTip
da	PICC perifert indført centralt katetersæt med Blue FlexTip kateter
nl	PICC Perifeer in te brengen centrale katheterset met Blue FlexTip katheter
fr	Jeu de cathéter central inséré périphériquement CCIP avec cathéter Blue FlexTip
de	PEZK-Besteck für peripher eingeführten, zentralen Katheter mit Blue-FlexTip-Katheter
el	Σετ περιφερικά εισαγόμενου κεντρικού καθετήρα PICC με καθετήρα Blue FlexTip
hu	PICC periferiálisán bevezethető, centrális katéterező felszerelés Blue FlexTip katéterrel
it	Set per catetere centrale a inserimento periferico con catetere Blue FlexTip
pl	Wprowadzany obwodowo zestaw typu PICC do cewnikowania naczyń centralnych z cewnikiem z końcówką Blue FlexTip
pt	Conjunto de Cateter Central de Introdução Periférica PICC com Cateter Blue FlexTip
ru	PICC набор центрального катетера для периферийного введения с катетером Blue FlexTip
sk	Sada periférne zavádzaného centrálného katétra (PICC) s katétrom Blue FlexTip
es	Conjunto de catéter de vena central de inserción periférica (PICC) con catéter Blue FlexTip
sv	Set för PICC perifert införd central kateter med Blue FlexTip-kateter

Diameter: 16 Ga., 1.7 mm	Priming Volume* (mL)	Flow Rate† (mL/hr)
Lumen		
Distal (16 Ga.)	0.64	1850

* Priming volumes are approximate and are done without accessories.

† Flow rates are with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IS-01651-130D (09/19)



PICC Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

REF PS-01651

16 Ga. 50 cm catheter length

ARROW



PICC Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

REF PS-01651

16 Ga. 50 cm catheter length

ARROW