

REF PS-01652

ARROW



PICC Two-Lumen Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

- cs** Dvoulumenový set s periferně zaváděným centrálním katetrem (PICC) s katetrem Blue FlexTip
- da** PICC perifert indført centralt katetersæt med dobbeltlumen og med Blue FlexTip kateter
- nl** PICC dubbellumen perifeer in te brengen centrale katheterset met Blue FlexTip katheter
- fr** Jeu de cathéter central inséré périphériquement CCIP à deux lumières avec cathéter Blue FlexTip
- de** PEZK-Besteck für zweilumigen, peripher eingeführten, zentralen Katheter mit Blue-FlexTip-Katheter
- el** Σετ περιφερικά εισαγόμενου κεντρικού καθετήρα δύο αυλών PICC με καθετήρα Blue FlexTip
- hu** PICC perifériálisan bevezethető, kétlumenű centrális katéterező felszerelés Blue FlexTip katéterrel
- it** Set per catetere centrale a inserimento periferico a due lumi con catetere Blue FlexTip
- pl** Wprowadzany obwodowo dwuświatłowy zestaw typu PICC do cewnikowania naczyń centralnych z cewnikiem z końcówką Blue FlexTip
- pt** Conjunto de Cateter Central de Introdução Periférica PICC com Lúmen Duplo e Blue FlexTip
- ru** PICC набор центрального двухпросветного катетера для периферийного введения с катетером Blue FlexTip
- sk** Sada dvojlúmenového periférne zavádzaného centrálneho katétra (PICC) s katétrom Blue FlexTip
- es** Conjunto de catéter de vena central de inserción periférica (PICC), de dos luces, con catéter Blue FlexTip
- sv** Set för PICC tvåkanalig perifert införd central kateter med Blue FlexTip-kateter

## Contents:

- 1: Two-Lumen Indwelling Catheter: 16 Ga. x 50 cm Radiopaque Polyurethane with Blue FlexTip®, Contamination Guard
- 1: Peel-Away Catheter: 14 Ga. x 2-1/2" (6.35 cm) Radiopaque over 15 Ga. TW Introducer Needle and 3 mL Luer-Slip Syringe
- 1: Skin Protectant Prep Pad
- 2: Dust Cap: Non-Vented
- 2: Clamp: Extension Line
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Dressing: STATLOCK®<sup>1</sup> Catheter Stabilization Device

IS-01652-130D (11/17)

<sup>1</sup>A registered trademark of C. R. Bard, Inc

All components are CE 0086 unless otherwise noted.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.40	575	3000
Proximal (20 Ga.)	0.30	160	1200

\* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates are determined at maximum pump pressure of 1.0 psig.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex  
Fluid path components are non-pyrogenic.



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

REF PS-01652

ARROW



PICC Two-Lumen Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

- cs** Dvoulumenový set s periferně zaváděným centrálním katetrem (PICC) s katetrem Blue FlexTip
- da** PICC perifert indført centralt katetersæt med dobbeltlumen og med Blue FlexTip kateter
- nl** PICC dubbellumen perifeer in te brengen centrale katheterset met Blue FlexTip katheter
- fr** Jeu de cathéter central inséré périphériquement CCIP à deux lumières avec cathéter Blue FlexTip
- de** PEZK-Besteck für zweilumigen, peripher eingeführten, zentralen Katheter mit Blue-FlexTip-Katheter
- el** Σετ περιφερικά εισαγόμενου κεντρικού καθετήρα δύο αυλών PICC με καθετήρα Blue FlexTip
- hu** PICC perifériálisan bevezethető, kétlumenű centrális katéterező felszerelés Blue FlexTip katéterrel
- it** Set per catetere centrale a inserimento periferico a due lumi con catetere Blue FlexTip
- pl** Wprowadzany obwodowo dwuświatłowy zestaw typu PICC do cewnikowania naczyń centralnych z cewnikiem z końcówką Blue FlexTip
- pt** Conjunto de Cateter Central de Introdução Periférica PICC com Lúmen Duplo e Blue FlexTip
- ru** PICC набор центрального двухпросветного катетера для периферийного введения с катетером Blue FlexTip
- sk** Sada dvojlúmenového periférne zavádzaného centrálneho katétra (PICC) s katétrom Blue FlexTip
- es** Conjunto de catéter de vena central de inserción periférica (PICC), de dos luces, con catéter Blue FlexTip
- sv** Set för PICC tvåkanalig perifert införd central kateter med Blue FlexTip-kateter

## Contents:

- 1: Two-Lumen Indwelling Catheter: 16 Ga. x 50 cm Radiopaque Polyurethane with Blue FlexTip®, Contamination Guard
- 1: Peel-Away Catheter: 14 Ga. x 2-1/2" (6.35 cm) Radiopaque over 15 Ga. TW Introducer Needle and 3 mL Luer-Slip Syringe
- 1: Skin Protectant Prep Pad
- 2: Dust Cap: Non-Vented
- 2: Clamp: Extension Line
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Dressing: STATLOCK®<sup>1</sup> Catheter Stabilization Device

IS-01652-130D (11/17)

<sup>1</sup>A registered trademark of C. R. Bard, Inc

All components are CE 0086 unless otherwise noted.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.40	575	3000
Proximal (20 Ga.)	0.30	160	1200

\* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates are determined at maximum pump pressure of 1.0 psig.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex  
Fluid path components are non-pyrogenic.



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com



PICC Two-Lumen Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

REF PS-01652 16 Ga. | 2 Lumen | 50 cm catheter length | ARROW



PICC Two-Lumen Peripherally Inserted Central Catheter Set with Blue FlexTip® Catheter

REF PS-01652 16 Ga. | 2 Lumen | 50 cm catheter length | ARROW