



REF CS-12402

4 Fr. | 2 Lumen | 5 cm catheter length | .018 inch dia. spring-wire guide

# Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**Contents:**

- 1: Two-Lumen Indwelling Catheter: 4 Fr. x 5 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.


Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65  
Warning: Cancer and Reproductive Harm [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central pédiatrique à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter für Kinder mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale a due lumi pediatrico con catetere Blue FlexTip
- pl** Dwukanałowy zestaw pediatryczny do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ru** Педиатрический двухпросветный комплект для катетеризации центральных вен с катетером с кончиком Blue FlexTip
- sl** Pediatrični komplet za centralno vensko katetrizacijo z dvema svetlinama s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces para uso pediátrico con catéter Blue FlexTip
- sv** Set för pediatrik tvåkanalig central venkatetrisering med Blue FlexTip-kateter
- tr** Pediyatrik İki Lümenli Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (22 Ga.)	0.24	930
Proximal (22 Ga.)	0.29	920

\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



CE 0086  
Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex, Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

IS-12402-161A Rev. 1 (08/18)



REF CS-12402

4 Fr. | 2 Lumen | 5 cm catheter length | .018 inch dia. spring-wire guide

Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

