

# ARROW®



## Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-12553-J** | **5.5** Fr. | **3** Lumen | **5** cm catheter length | **.018** inch dia. spring-wire guide

### Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 5 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 3: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrering med Blue FlexTip-kateter




Not made with natural rubber latex

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (20 Ga.)	0.27	2000
Medial (22 Ga.)	0.29	1400
Proximal (22 Ga.)	0.33	1400

\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

**CE** REP 0086  
Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Teleflex®**  
www.teleflex.com  
Manufacturer:  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

IS-12553-160C (01/18)

DuPont™ Tyvek®  
production



## Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-12553-J**

**5.5** Fr. | **3** Lumen | **5** cm catheter length | **.018** inch dia. spring-wire guide

# ARROW®