



Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-12702-E**

7 Fr. | **2** Lumen | **16** cm catheter length | **.025** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- sv** Set för tvåkanalig central venkatetrisering med Blue FlexTip-kateter


All components are CE 0086 unless otherwise noted.



Rx only.

Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (18 Ga.)	0.33	1500
Proximal (14 Ga.)	0.52	5500

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA



IS-12702-130B (01/20)



Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-12702-E**

7 Fr. | **2** Lumen | **16** cm catheter length | **.025** inch dia. spring-wire guide