



Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **ES-14402**

4 Fr. | **2** Lumen | **13** cm catheter length | **.018** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 4 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 22 Ga. x 1-1/16" (2.67 cm) TW
- 1: Syringe: 3 mL Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

All components are CE 0086 unless otherwise noted.




Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central pédiatrique à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter für Kinder mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale a due lumi pediatrico con catetere Blue FlexTip
- pl** Dwukanałowy zestaw pediatryczny do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ru** Педиатрический двухпросветный комплект для катетеризации центральных вен с катетером с кончиком Blue FlexTip
- sl** Pediatrični komplet za centralno venosko katetrizacijo z dvema svetlinama s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces para uso pediátrico con catéter Blue FlexTip
- sv** Set för pediatrik tvåkanalig central venkatetrering med Blue FlexTip-kateter
- tr** Blue FlexTip Kateterli Pediyatrik İki Lümenli Santral Venöz Kateterizasyon Seti



Lumen	Priming Volume *(mL)	Flow Rate† (mL/hr)
Distal (22 Ga.)	0.25	550
Proximal (22 Ga.)	0.31	425

* Priming volumes are approximate and are done without accessories.
 † Based on 100 cm head height.



Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IS-14402-162G (08/19)



Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **ES-14402**

4 Fr. | **2** Lumen | **13** cm catheter length | **.018** inch dia. spring-wire guide