

# ARROW®



## Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-14502**

**5** Fr. | **2** Lumen | **13** cm catheter length | **.021** inch dia. spring-wire guide

**Contents:**

- 1: Two-Lumen Indwelling Catheter: 5 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .021" (0.53 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW
- 1: Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central pédiatrique à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter für Kinder mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale a due lumi pediatrico con catetere Blue FlexTip
- pl** Dwukanałowy zestaw pediatryczny do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ru** Педиатрический двухпросветный комплект для катеризации центральных вен с катетером с кончиком Blue FlexTip
- sl** Pediatrični komplet za centralno vensko kateterizacijo z dvema svetlinama s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces para uso pediátrico con catéter Blue FlexTip
- sv** Set för pediatriisk tvåkanalig central venkateterisering med Blue FlexTip-kateter
- tr** Blue FlexTip Kateterli Pediyatrik İki Lümenli Santral Venöz Kateterizasyon Seti




Not made with natural rubber latex

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (18 Ga.)	0.30	1500
Proximal (20 Ga.)	0.30	760

\* Priming volumes are approximate and are done without accessories.  
 † Based on 100 cm head height.

DuPont™ Tyvek® packaging

**Manufacturer:**  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**CE** **REP** Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
 0086

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IS-14502-162A (03/18)



## Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-14502**

**5** Fr. | **2** Lumen | **13** cm catheter length | **.021** inch dia. spring-wire guide

# ARROW®