



Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for Femoral Vein Insertion

REF **CS-14553** **5.5** Fr. | **3** Lumen | **30** cm catheter length | **.018** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 30 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.




Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (20 Ga.)	0.38	800
Medial (22 Ga.)	0.35	300
Proximal (22 Ga.)	0.39	300

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- fr** Jeu de cathétérisme veineux central pédiatrique à lumières multiples avec cathéter Blue FlexTip pour insertion dans la veine fémorale
- de** Besteck für viellumigen, zentralen Venenkatheter für Kinder mit Blue-FlexTip-Katheter zur Einführung in die Femoralvene
- it** Set per cateterismo della vena centrale multilume pediatrico con catetere Blue FlexTip per introduzione nella vena femorale
- pl** Wielokanałowy zestaw pediatryczny do cewnikowania żył centralnych z cewnikiem Blue FlexTip, do wprowadzania przez żyłę udową
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, para introdução na Veia Femoral
- ru** Педиатрический многопросветный комплект для катетеризации центральных вен с катетером с кончиком Blue FlexTip для введения в бедренную вену
- sl** Pediatrični komplet za centralno venosko katetizacijo z več svetlinami s katetrom Blue FlexTip za vstavitve v stegensko veno
- es** Conjunto de cateterización venosa central de luz múltiple para uso pediátrico con catéter Blue FlexTip para inserción en la vena femoral
- sv** Set för pediatriisk flerkanalig central venkatetrering med Blue FlexTip-kateter för införing i lårbensven
- tr** Femoral Ven İnseriyonu için Blue FlexTip Kateterli Pediyatrik Çoklu Lümenli Santral Venöz Kateterizasyon Seti



Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IS-14553-162A (03/18)

DuPont™ Tyvek®
production



Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for Femoral Vein Insertion

REF **CS-14553**

5.5 Fr. | **3** Lumen | **30** cm catheter length | **.018** inch dia. spring-wire guide