

ARROW®



Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-14703-E** | 7^{Fr.} | 3^{Lumen} | 30^{cm} catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 30 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**


California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- pl** Wielokanałowy zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ru** Комплект для катетеризации главных вен с многопросветным катетером Blue FlexTip
- sl** Komplet za centralno vensko katetrizacijo z več svetlinami s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip-kateter
- tr** Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



Not made with natural rubber latex



* Priming volumes are approximate and are done without accessories.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IS-14703-200A (05/18)



Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-14703-E**

7^{Fr.} | 3^{Lumen} | 30^{cm} catheter length | **.032** inch dia. spring-wire guide

ARROW®