



## Large-Bore Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

REF **CS-15123-F**      **12** Fr. | **3** Lumen | **20** cm catheter length | **.035** inch dia. spring-wire guide

### Contents:

- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 3: Dust Cap: Non-Vented
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11


All components are CE 0086 unless otherwise noted.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)
Distal (16 Ga.)	0.40	80
Medial (12 Ga.)	1.60	285
Proximal (12 Ga.)	1.50	305

\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- fr** Jeu de cathétérisme veineux central de grand calibre à lumières multiples avec cathéter Blue FlexTip pour les perfusions de grand volume
- de** Besteck für viellumigen, zentralen, großkalibrigen Venenkatheter mit Blue-FlexTip-Katheter für Infusionen großer Flüssigkeitsmengen
- it** Set per cateterismo della vena centrale multilume di diametro ampio con catetere Blue FlexTip per volumi di infusione elevati
- pl** Wielokanałowy zestaw o dużej średnicy do cewnikowania żył centralnych z cewnikiem Blue FlexTip do przetwarzania dużych objętości płynów infuzyjnych
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Grande Calibre, com Cateter Blue FlexTip, para Infusões de Grandes Volumes
- ru** Комплект для катетеризации главных вен с многопросветным катетером большого внутреннего диаметра с Blue FlexTip для вливаний большого объема
- sl** Komplet za centralno vensko katetrizacijo z velikimi luknjami in več svetlinami s katetrom Blue FlexTip za infundiranje velike prostornine
- es** Conjunto de cateterización venosa central de luz múltiple y gran diámetro con catéter Blue FlexTip para infusiones de gran volumen
- sv** Set för storhålig flerkanalig central venkatetrering med Blue FlexTip-kateter för infusioner med hög volym
- tr** Büyük Açıklıklı Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile, Yüksek Hacimli İnfüzyonlar İçin



STERILE EO



Not made with natural rubber latex

**Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IS-15123-135A (05/18)



## Large-Bore Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

REF **CS-15123-F**

**12** Fr. | **3** Lumen | **20** cm catheter length | **.035** inch dia. spring-wire guide