

ARROW®



Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-15553-E** | **5.5** Fr. | **3** Lumen | **8** cm catheter length | **.018** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 8 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central pédiatrique à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter für Kinder mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale multilume pediatrico con catetere Blue FlexTip
- pl** Wielokanałowy zestaw pediatryczny do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ru** Педиатрический многопросветный комплект для катетеризации центральных вен с катетером с кончиком Blue FlexTip
- sl** Pediatrični komplet za centralno vensko katetризациjo z več svetlinami s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple para uso pediátrico con catéter Blue FlexTip
- sv** Set för pediatrik flerkanalig central venkatetrisering med Blue FlexTip-kateter
- tr** Pediyatrik Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



STERILE EO




Not made with natural rubber latex

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (20 Ga.)	0.28	1300
Medial (22 Ga.)	0.30	900
Proximal (22 Ga.)	0.32	1050

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

CE **REP** Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 0086

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA
Teleflex®
www.teleflex.com

IS-15553-160D (01/18)

DuPont™ Tyvek®
 packaging



Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-15553-E** | **5.5** Fr. | **3** Lumen | **8** cm catheter length | **.018** inch dia. spring-wire guide

ARROW®