



## Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-15703-E** | **7** Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

### Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm**  
**www.P65Warnings.ca.gov.**

Fluid path components are non-pyrogenic.


- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- pl** Wielokanałowy zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ru** Комплект для катетеризации главных вен с многопросветным катетером Blue FlexTip
- sl** Komplet za centralno vensko katetrizacijo z več svetlinami s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrering med Blue FlexTip-kateter
- tr** Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



STERILE EO



Not made with natural rubber latex



\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

**CE** REP 0086 Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IS-15703-137A (05/18)



## Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-15703-E**

**7** Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide