

ARROW®



Arrow-Howes™ Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-15854-E** | **8.5** Fr. | **4** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central à quatre lumières Arrow-Howes avec cathéter Blue FlexTip
- de** Arrow-Howes-Besteck für vierlumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale a quattro lumi Arrow-Howes con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central Arrow-Howes com Lúmen Quádruplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de cuatro luces Arrow-Howes con catéter Blue FlexTip
- sv** Arrow-Howes set för fyrkanalig central venkatetrisering med Blue FlexTip-kateter




Not made with natural rubber latex

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.43	2700
Medial 1 (14 Ga.)	0.54	5100
Medial 2 (18 Ga.)	0.36	1400
Proximal (18 Ga.)	0.41	1600

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

CE **REP** Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 0086

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IS-15854-130B (01/18)

DuPont™ Tyvek®
 packaging



Arrow-Howes™ Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-15854-E** | **8.5** Fr. | **4** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

ARROW®