



REF CS-16553-E

5.5 Fr. | 3 Lumen | 13 cm catheter length | .018 inch dia. spring-wire guide

Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp


All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central pédiatrique à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter für Kinder mit Blue-Flex-Tip-Katheter
- it** Set per cateterismo della vena centrale multilume pediatrico con catetere Blue FlexTip
- pl** Wielokanałowy zestaw pediatriczny do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ru** Педиатрический многопросветный комплект для катетеризации центральных вен с катетером с кончиком Blue FlexTip
- sl** Pediatrični komplet za centralno vensko katetrizacijo z več svetlinami s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple para uso pediátrico con catéter Blue FlexTip
- sv** Set för pediatrik flerkanalig central venkatetrering med Blue FlexTip-kateter
- tr** Pediyatrik Çoklu Lümenli Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (20 Ga.)	0.30	1200
Medial (22 Ga.)	0.31	800
Proximal (22 Ga.)	0.34	900

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

CE 0086
Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

IS-16553-161A (01/18)



REF CS-16553-E 5.5 Fr. | 3 Lumen | 13 cm catheter length | .018 inch dia. spring-wire guide

Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

