

Two-Lumen Central Venous Catheterization Set

with Blue FlexTip® Catheter




ARROW
INTERNATIONAL

Arrow Select Kits

REF **CV-17702-S**

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps, Injection Site Cap
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: Pressure Transduction Probe
- 1: Safety Introducer Needle: 18 Ga. x 2-3/4" (7 cm) Echogenic
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator
- 1: Dust Cap: Non-Vented
- 1: Clamp: Catheter
- 1: Fastener: Catheter Clamp



Lumen	Priming Vol.* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.57	5000
Proximal (18 Ga.)	0.40	1500

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road, Reading, PA 19605 USA
www.arrowintl.com

IS-17702-109A (6/12)

Two-Lumen Central Venous Catheterization Set

REF **CV-17702-S**

7 Fr.

2 Lumen

20 cm catheter length

.032 in. dia. spring-wire guide

Arrow Select Kits



ARROW
INTERNATIONAL