


Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	1.05	2500
Proximal (18 Ga.)	0.71	590

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 24" (60 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps, Injection Site Caps
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 39-1/2" (100 cm) PTFE-Coated (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator
- 1: Clamp: Catheter
- 1: Fastener: Catheter Clamp

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use — do not sterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

Arrow International, Inc.
2400 Bernville Road, Reading, PA 19605 USA

TELEFLEX Medical
IDA Business and Technology Park, Athlone, Ireland



REF CS-17752

Two-Lumen Antecubital Central Venous Catheterization Set with Blue FlexTip® Catheter



- FR** Jeu de cathétérisme veineux central antébrachial à deux lumières avec cathéter Blue FlexTip
- DE** Besteck für zweilumigen, zentralen Venenkatheter zur Einführung in der Ellenbeuge mit Blue-FlexTip-Katheter
- IT** Set per cateterismo della vena centrale cubitale anteriore a due lumi con catetere Blue FlexTip
- PL** Dwukanałowy łokciowy zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central Antecubital com Lúmen Duplo e Cateter Blue FlexTip
- RU** Комплект для катетеризации главных вен доступом спереди от локтевого сустава с двухпросветным катетером Blue FlexTip
- SL** Komplet za antekubitalno centralno vensko katetrizacijo z dvema svetlinama s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central antecubital de dos luces con catéter Blue FlexTip
- SV** Set för tvåkanalig antekubital central venkatetrering med Blue FlexTip-kateter
- TR** İki Lümenli Antekübital Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile

15-17752-1308 (12/11)



REF
CS-17752

7 Fr.

2 Lumen

60 cm catheter length

.035 inch dia. spring-wire guide

Two-Lumen Antecubital Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open