

# ARROW®



## Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

REF **CS-24402**

**4** Fr. | **2** Lumen | **13** cm catheter length | **.018** inch dia. spring-wire guide

**Contents:**

- 1: Two-Lumen Indwelling Catheter: 4 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

<sup>1</sup>Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.




Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.

**Indications For Use:** • Provides central venous access. • Provides protection against catheter related infection. • 4 Fr. antimicrobial catheter should be used in pediatric patients weighing 2.0 kg or greater.

- fr** Jeu de cathétérisme veineux central pédiatrique à deux lumières avec cathéter Blue FlexTip ARROWgard Blue
- de** Besteck für zweilumigen, zentralen Venenkatheter für Kinder mit Blue-FlexTip-ARROWgard-Blue-Katheter
- it** Set per cateterismo della vena centrale a due lumi pediatrico con catetere Blue FlexTip ARROWgard Blue
- pl** Dwukanałowy zestaw pediatryczny do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWgard Blue
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWgard Blue
- ru** Педиатрический комплект для катетеризации главных вен с двухпросветным катетером Blue FlexTip ARROWgard Blue
- sl** Pediatrični komplet za centralno vensko katetizacijo z dvema svetlinama s katetrom Blue FlexTip ARROWgard Blue
- es** Conjunto de cateterización venosa central de dos luces para uso pediátrico con catéter Blue FlexTip ARROWgard Blue
- sv** Set för pediatrik tvåkanalig central venkatetrisering med Blue FlexTip ARROWgard Blue-kateter
- tr** Pediyatrik İki Lümenli Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWgard Blue Kateter ile

**Contraindications:** The ARROWgard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (22 Ga.)	0.27	620
Proximal (22 Ga.)	0.32	460

\* Priming volumes are approximate and are done without accessories.  
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

Manufacturer:  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA



IS-24402-200A (05/18)



## Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

REF **CS-24402**

**4** Fr. | **2** Lumen | **13** cm catheter length | **.018** inch dia. spring-wire guide

# ARROW®