



## Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter

REF CS-25703-E

7<sup>Fr.</sup> | 3 Lumen | 20 cm catheter length | .032 inch dia. spring-wire guide

### Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup>Licensed under US Patent No. 6,706,024.


All components are CE 0086 unless otherwise noted.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.

**Contraindications: The ARROWgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).**

- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWgard Blue
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-ARROWgard-Blue-Katheter
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWgard Blue
- pl** Wielokanałowy zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWgard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWgard Blue
- ru** Комплект для катетеризации главных вен с многопросветным катетером Blue FlexTip ARROWgard Blue
- sl** Komplet za centralno vensko katetrizacijo z več svetlinami s katetrom Blue FlexTip ARROWgard Blue
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWgard Blue
- sv** Set för flerkanalig central venkatetrering med Blue FlexTip ARROWgard Blue-kateter
- tr** Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWgard Blue Kateter ile



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.40	2720
Medial (18 Ga.)	0.40	1260
Proximal (18 Ga.)	0.40	1390

\* Priming volumes are approximate and are done without accessories.  
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



CE 0086 Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

Manufacturer:  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**Teleflex®**  
 www.teleflex.com

IS-25703-131A (06/18)



Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter  
 REF CS-25703-E

7<sup>Fr.</sup> | 3 Lumen | 20 cm catheter length | .032 inch dia. spring-wire guide