



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

REF LU-09803-S

8.5 Fr. | 10 cm sheath length | .025 inch dia. spring-wire guide

Contents:

- | | |
|--|--|
| <ul style="list-style-type: none"> 1: Sheath: 8.5 Fr. x 10 cm Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port, Tissue Dilator and Bonded High-Flow 3-way Stopcock 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle 1: Catheter: 16 Ga. x 1-7/8" (4.8 cm) Radiopaque over 19 Ga. TW Introducer Needle 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW 1: Local Anesthetic Component 2: Syringe: 5 mL Luer-Slip 1: 5 mL 1% Lidocaine HCl Solution | <ul style="list-style-type: none"> 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution) 1: Obturator: 8 Fr. 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape 1: SharpsAway® Disposal Cup 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration 1: Scalpel: #11 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) 1: Suture: 3-0 Braided Silk with Straight Needle 1: Tubing[‡]: 24" (61 cm) Extension |
|--|--|

¹A registered trademark of CareFusion or one of its subsidiaries.

							This product contains dry natural rubber.		‡ Precaution: Contains Phthalate: DEHP		Store between 20 - 25°C (68 - 77°F).	

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

DuPont™ Tyvek®
packaging

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IIU-09803-102M (12/17)



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

REF LU-09803-S

8.5 Fr. | 10 cm sheath length | .025 inch dia. spring-wire guide