



You-Bend™ Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Arrowg+ard Blue® Catheter for High Volume Infusions

REF **CU-23122-F**

12 Fr. | **2** Lumen | **13** cm catheter length | **.035** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter†: 12 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11

¹Licensed under US Patent No. 6,706,024.


All components are CE 0086 unless otherwise noted.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels. The Arrowg+ard Blue antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).


- fr** Jeu de cathétérisme pour hémodialyse à deux lumières You-Bend avec cathéter Blue FlexTip Arrowg+ard Blue pour les perfusions de grand volume
- de** Besteck für zweilumigen You-Bend-Hämodialysekatheter mit Blue-FlexTip-Arrowg+ard-Blue-Katheter für Infusionen großer Flüssigkeitsmengen
- it** Set per cateterismo in emodialisi a due lumi You-Bend con catetere Blue FlexTip Arrowg+ard Blue per volumi di infusione elevati
- pl** Dwukanałowy zestaw You-Bend do hemodializy z cewnikiem Blue FlexTip Arrowg+ard Blue do przetaczania dużych objętości płynów infuzyjnych
- pt** Conjunto de Cateterização para Hemodiálise You-Bend com Lúmen Duplo e Cateter Blue FlexTip, Arrowg+ard Blue, para Infusões de Grandes Volumes
- ru** Двухпросветный комплект для гемодиализной катетеризации You-Bend с катетером Blue FlexTip Arrowg+ard Blue для вливания больших объемов
- sl** Komplet za hemodializno katetризacijo You-Bend z dvema svetlinama s katetrom Blue FlexTip Arrowg+ard Blue za infundiranje velike prostornine
- es** Conjunto de cateterización de dos luces para hemodiálisis You-Bend con catéter Blue FlexTip Arrowg+ard Blue para infusiones de gran volumen
- sv** You-Bend set för tvåkanalig hemodialyskatetrisering med Blue FlexTip Arrowg+ard Blue-kateter för infusioner med hög volym
- tr** You-Bend İki Lümenli Hemodiyaliz Kateterizasyon Seti, Blue FlexTip Arrowg+ard Blue Kateterli, Yüksek Hacimli İnfüzyonlar için



Not made with natural rubber latex

‡ MR Conditional

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (12 Ga.)	1.50	230 400	100 200
Proximal (12 Ga.)	1.40	230 400	100 200

* Priming volumes are approximate and are done without accessories.
† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA



Teleflex®
www.teleflex.com

IU-23122-135A (2020-06)



You-Bend™ Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Arrowg+ard Blue® Catheter for High Volume Infusions

REF **CU-23122-F**

12 Fr. | **2** Lumen | **13** cm catheter length | **.035** inch dia. spring-wire guide