

Contents:

1:	Catheter: 20 Ga. x 4.45 cm Radiopaque over 22 Ga. TW Introducer Needle
1:	Wing Clip: Suture

Contents:

1:	Catheter: 20 Ga. x 4.45 cm Radiopaque over 22 Ga. TW Introducer Needle
1:	Wing Clip: Suture

Contents:

1: Catheter: 20 Ga. x 4.45 cm Radiopaque
over 22 Ga. TW Introducer Needle

1: Wing Clip: Suture

Contents:

1:	Catheter: 20 Ga. x 4.45 cm Radiopaque over 22 Ga. TW Introducer Needle
1:	Wing Clip: Suture

Contents:

1: Catheter: 20 Ga. x 4.45 cm Radiopaque
over 22 Ga. TW Introducer Needle
1: Wing Clip: Suture

1:	Catheter: 20 Ga. x 4.45 cm Radiopaque
1:	over 22 Ga. TW Introducer Needle
1:	Wing Clip: Suture

Contents:

1. Catheter: 20 Ga. x 4.5 cm Radiopaque
over 22 Ga. TW Introducer Needle

1. Wing Clip: Suture

Contents:

1. Catheter: 20 Ga. x 4.5 cm Radiopaque over 22 Ga. TW Introducer Needle

1. Wing Clip: Suture