

Catheter Diameter: 16 Ga., 1.7 mm

Priming Volume: 0.43 mL

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

Flow Rate: 2900 mL/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Contents:

- 1: . . . Indwelling Catheter: 16 Ga. x 8" (20 cm)
Radiopaque Polyurethane
with Blue FlexTip®, Integral Extension Line,
Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide, Marked:
.032" (0.81 mm) dia. x 23-5/8" (60 cm)
(Straight Soft Tip on One End—
"J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle:
18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator
- 1: Clamp: Catheter
- 1: Fastener: Catheter Clamp

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use — do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

 Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Athlone, Ireland



REF CV-04301

**Central Venous
Catheterization Set
with Blue FlexTip® Catheter**

ARROW®

- FR** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- DE** Besteck für zentralen Venenkatheter mit Blue-FlexTip-Katheter
- IT** Set per cateterismo della vena centrale con catetere Blue FlexTip
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- RU** Комплект для катетеризации главных вен с катетером Blue FlexTip
- SL** Komplet za centralno venosko katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central con catéter Blue FlexTip
- SV** Set för central venkatetrisering med Blue FlexTip-kateter
- TR** Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile

IN-04301-1308 (9/11)

ARROW®
INTERNATIONAL

REF
CV-04301

16^{Ga.} **20** cm catheter length **.032** inch dia. spring-wire guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open