

Catheter Diameter: 16 Ga., 1.7 mm

Priming Volume: 0.40 mL

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

Flow Rate: 3700 mL/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.


**Contents:**

- 1: ..... Indwelling Catheter:  
16 Ga. x 6" (16 cm) Radiopaque Polyurethane  
with Blue FlexTip®, Integral Extension Line,  
Extension Line Clamp, Integral Suture Wing
- 1: ..... Spring-Wire Guide, Marked:  
.032" (0.81 mm) dia. x 17-3/4" (45 cm)  
(Straight Soft Tip on One End –  
"J" Tip on Other) with Arrow Advancer™
- 1: ..... Introducer Needle:  
18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: ..... Syringe: 5 mL Luer-Slip
- 1: ..... Tissue Dilator
- 1: ..... Clamp: Catheter
- 1: ..... Fastener: Catheter Clamp

**Rx only.**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

 Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

 Teleflex Medical  
IDA Business and Technology Park  
Athlone, Ireland



**REF CV-04306**

**Central Venous  
Catheterization Set  
with Blue FlexTip® Catheter**

**ARROW®**

- FR** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- DE** Besteck für zentralen Venenkatheter mit Blue-FlexTip-Katheter
- IT** Set per cateterismo della vena centrale con catetere Blue FlexTip
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- RU** Комплект для катеризации главных вен с катетером Blue FlexTip
- SL** Komplet za centralno venosko katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central con catéter Blue FlexTip
- SV** Set för central venkatetrisering med Blue FlexTip-kateter
- TR** Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile

IV04306-1308 (3/12)

**ARROW®**  
INTERNATIONAL

REF  
**CV-04306**

**16**<sup>Ga.</sup> **16**<sup>cm</sup> catheter length **.032**<sup>inch dia.</sup> spring-wire guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to  
Open