



Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

REF CV-04701-BZ

pt Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
Manter afastado da luz solar. Manter seco.

es Conjunto de cateterización venosa central con catéter Blue FlexTip
Mantener alejado de la luz del sol. Mantener en un lugar seco.

Contents:

- 1: Indwelling Catheter: 14 Ga. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Side Holes, Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.



Caution



Consult instructions for use



Do not reuse



Do not resterilize

STERILE EO

Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

Manufacturer
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

CE **REF** Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

Catheter Diameter: 14 Ga., 2.2 mm

Priming Volume: 0.60 mL

Priming volumes are approximate and are done without accessories.

Flow Rate: 6880 mL/hr

Based on 100 cm head height.

IV-04701-131A (3/17)



Central Venous Catheterization Set with Blue FlexTip® Catheter

REF CV-04701-BZ

14 Ga. | 20 cm catheter length | .032 inch dia. spring-wire guide

ARROW

