

ARROW®



Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CV-15553-BZ** 5.5 Fr. | 3 Lumen | 8 cm catheter length | **.018** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 8 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

pt Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip

es Conjunto de cateterización venosa central de luz múltiple para uso pediátrico con catéter Blue FlexTip


All components are CE 0086 unless otherwise noted.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (20 Ga.)	0.28	1300
Medial (22 Ga.)	0.30	900
Proximal (22 Ga.)	0.32	1050

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

Catheter Cross Section

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
CE 0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IV-15553-161B (11/19)



Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CV-15553-BZ**

5.5 Fr. | 3 Lumen | 8 cm catheter length | **.018** inch dia. spring-wire guide

ARROW®