

ARROWgard Blue PLUS⁺ Multi-Lumen CVC Kit with Blue FlexTip[®] Catheter



Arrow Select Kits

ASK-42703-RV1

- > **AGB⁺ Catheter**
- > **Sharps Safety Features**
- > **Maximal Barrier Precautions**

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip[®], ARROWgard Blue PLUS⁺ Antimicrobial Surface Treatment¹, Extension Line Clamps**
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer[™]**
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle**
- 1: Introducer Needle: 18 Ga. x 1-1/2" (3.81 cm) XTW**
- 1: Pressure Transduction Probe**
- 1: Syringe: 10 mL Luer-Lock**
- 1: Injection Needle: SafetyGlide^{®2} 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe**
- 1: Injection Needle: SafetyGlide^{®2} 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe**
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe**
- 1: Tissue Dilator**
- 1: Suture: 3-0 Silk with Curved Needle**
- 1: SecondSite[™] Adjustable Hub: Catheter Clamp and Fastener**
- 1: SharpsAway II[™] Locking Disposal Cup**
- 2: CSR Wraps**
- 1: Maximal Barrier Drape[™]: 54" x 96" with 4" fenestration**
- 1: Safety Scalpel: #11**
- 4: Staple Anchoring Devices^{3,4}**
- 1: Needle Holder**
- 1: Forceps**
- 2: Gauze Pads: 2" x 2"**
- 5: Gauze Pads: 4" x 4"**
- 3: Prep Sponge Swabs**
- 1: Surgical Apparel: Bouffant Cap**
- 1: Surgical Apparel: Impervious Gown**
- 1: Surgical Apparel: Mask with Eye Shield**
- 1: Tray: Prep**

¹ Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

² A registered trademark of Becton, Dickinson and Company.

³ Licensed under US Patent No. RE37,814.

Infection Protection Feature:

- > **AGB⁺ catheters** include internal lumen protection & external antimicrobial surfaces utilizing chlorhexidine, chlorhexidine acetate and silver sulfadiazine

Maximal Barrier Precautions:

- > Drape Sheet
- > Bouffant Surgical Cap
- > Impervious Surgical Gown
- > Surgical Mask with Eye Shield

Safety Features:

- > SafetyGlide[®] Protected Needles
- > SharpsAway II[™] Locking Disposal Cup
- > Safety Scalpel
- > Staple Anchoring Devices

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Contraindications: The ARROWgard Blue PLUS⁺ antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, chlorhexidine acetate, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).


Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication.

⁴MR Conditional.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).



Arrow International, Inc.
2400 Bernville Road, Reading, PA 19605 USA
www.arrowintl.com



Lumen	Priming Volume* (mL)	Flow Rate [†] (mL/hr)
Distal (16 Ga.)	0.38	3000
Medial (18 Ga.)	0.37	1400
Proximal (18 Ga.)	0.39	1600

*Priming volumes are approximate and are done without an injection cap.

[†] Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

ARROWgard Blue PLUS⁺ Multi-Lumen CVC Kit

ASK-42703-RV1

Arrow Select Kits



7 Fr.

3 Lumen

16 cm catheter length

.032 inch diameter spring-wire guide

- > **AGB⁺ Catheter**
- > **Sharps Safety Features**
- > **Maximal Barrier Precautions**

IV-42703-104A (1/14)