

ARROW®



Arrowg+ard Blue PLUS® Multi-Lumen CVC Set

REF **CV-45703-E**

7^{Fr.} | 3^{Lumen} | 20^{cm} catheter length | .032^{inch dia.} spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard Blue PLUS® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent Nos. 6,706,024, 6,872,195 and 7,329,412.


All components are CE 0086 unless otherwise noted.



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.44	3100
Medial (18 Ga.)	0.39	1500
Proximal (18 Ga.)	0.39	1600

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- cs** Vicolumenový CVC set Arrowg+ard Blue PLUS
- nl** Arrowg+ard Blue PLUS multilumen centraal-veneuzen katheterset
- fr** Jeu de cathéters veineux centraux à lumières multiples Arrowg+ard Blue PLUS
- de** Set mit mehrlumigen Arrowg+ard Blue PLUS ZVK
- it** Set CVC multilume Arrowg+ard Blue PLUS
- pl** Wielokanałowy zestaw do cewnikowania żył centralnych z cewnikiem Arrowg+ard Blue PLUS
- pt** Conjunto de CVC de Lúmen Múltiplo Arrowg+ard Blue PLUS
- ru** Комплект для катетеризации главных вен с многопросветным катетером Arrowg+ard Blue PLUS
- sl** Komplet za centralno vensko katetrizacijo z več svetlinami s katetrom Arrowg+ard Blue PLUS
- es** Conjunto de CVC de luz múltiple Arrowg+ard Blue PLUS
- sv** Arrowg+ard Blue PLUS flerlumen CVK-set
- tr** Multi-Lümen Santral Venöz Kateterizasyon Seti, Arrowg+ard Blue PLUS Kateter ile

EC REP
CE 0086
Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IV-45703-134C (09/19)



**Arrowg+ard Blue PLUS®
Multi-Lumen CVC Set**
REF **CV-45703-E**

7^{Fr.} | 3^{Lumen} | 20^{cm} catheter length | .032^{inch dia.} spring-wire guide

ARROW®