

ARROW®



Arrowgard Blue PLUS® Two-Lumen CVC Set

REF **CV-45802-E**

8 Fr. | **2** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 8 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowgard Blue Plus® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent Nos. 6,706,024, 6,872,195 and 7,329,412.


All components are CE 0086 unless otherwise noted.



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.

Contraindications: The Arrowgard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.75	5600
Proximal (14 Ga.)	0.80	5100

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- cs** Dvoulumenový CVC set Arrowgard Blue PLUS
- nl** Arrowgard Blue PLUS dubbellumen centraal-veneuze katheterset
- fr** Jeu de cathéters veineux centraux à deux lumières Arrowgard Blue PLUS
- de** Set mit zweilumigen Arrowgard Blue PLUS ZVK
- it** Set CVC bilume Arrowgard Blue PLUS
- pl** Dwukanałowy zestaw do cewnikowania żył centralnych z cewnikiem Arrowgard Blue PLUS
- pt** Conjunto de CVC de Lúmen Duplo Arrowgard Blue PLUS
- ru** Комплект для катетеризации главных вен с двухпросветным катетером Arrowgard Blue PLUS
- si** Komplet za centralno vensko katetrizacijo z dvema svetlinama s katetrom Arrowgard Blue PLUS
- es** Conjunto de CVC de doble luz Arrowgard Blue PLUS
- sv** Arrowgard Blue PLUS CVK-set med två lumen
- tr** İki-Lümenli Santral Venöz Kateterizasyon Seti, Arrowgard Blue PLUS Kateter ile

CE **REP** Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 0086

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
 www.teleflex.com

IV-45802-130C (09/19)



Arrowgard Blue PLUS® Two-Lumen CVC Set

REF **CV-45802-E**

8 Fr. | **2** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

ARROW®