

REF CV-50016

Central Venous Catheterization Set

ARROW®

Catheter Diameter: 16 Ga., 1.7 mm

Priming Volume: 0.22 mL

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

Flow Rate: 3400 mL/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Contents:

- 1: Indwelling Catheter: 16 Ga. x 8" (20 cm) Radiopaque Polyurethane with Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator
- 1: Clamp: Catheter
- 1: Fastener: Catheter Clamp

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use — do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

 Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Athlone, Ireland



- FR** Jeu de cathétérisme veineux central
- DE** Besteck für zentralen Venenkatheter
- IT** Set per cateterismo della vena centrale
- PL** Zestaw do cewnikowania żył centralnych
- PT** Conjunto de Cateterização Venosa Central
- RU** Комплект для катетеризации главных вен
- SL** Komplet za centralno venosko katetrizacijo
- ES** Conjunto de cateterización venosa central
- SV** Set för central venkatetrering
- TR** Santral Venöz Kateterizasyon Seti

IV-50016-130C (08/15)

ARROW®

INTERNATIONAL

REF CV-50016 16^{Ga.} 20^{cm} catheter length .032^{inch dia.} spring-wire guide
Central Venous Catheterization Set

Peel to Open