

ARROW®



Spring-Wire Guide .018" (0.46 mm) diameter x 51-3/16" (130 cm) Marked Wire

REF **PW-18130-P**

.018 inch dia. spring-wire guide | **130** cm spring-wire guide

Contents:

- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 51-3/16" (130 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

DuPont™ Tyvek®
packaging

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

 **California Prop. 65**

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



STERILE EO



Not made with natural rubber latex

 **Manufacturer:**
Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

IW-18130-101E, Rev. 01 (2020-12)



Spring-Wire Guide .018" (0.46 mm) diameter x 51-3/16" (130 cm) Marked Wire

REF **PW-18130-P**

.018 inch dia. spring-wire guide | **130** cm spring-wire guide

ARROW®